








Tube #	Salt	Tube #	Polymer	Tube #	pH \diamond	
1.	0.2 M Sodium fluoride	1.	20% w/v Polyethylene glycol 3,350	1.	7.3	F ⁻
2.	0.2 M Potassium fluoride	2.	20% w/v Polyethylene glycol 3,350	2.	7.3	Cl ⁻
3.	0.2 M Ammonium fluoride	3.	20% w/v Polyethylene glycol 3,350	3.	6.2	I ⁻
4.	0.2 M Lithium chloride	4.	20% w/v Polyethylene glycol 3,350	4.	6.8	Fluoride Chloride Iodide
5.	0.2 M Magnesium chloride hexahydrate	5.	20% w/v Polyethylene glycol 3,350	5.	5.9	
6.	0.2 M Sodium chloride	6.	20% w/v Polyethylene glycol 3,350	6.	6.9	
7.	0.2 M Calcium chloride dihydrate	7.	20% w/v Polyethylene glycol 3,350	7.	5.1	
8.	0.2 M Potassium chloride	8.	20% w/v Polyethylene glycol 3,350	8.	7.0	
9.	0.2 M Ammonium chloride	9.	20% w/v Polyethylene glycol 3,350	9.	6.3	
10.	0.2 M Sodium iodide	10.	20% w/v Polyethylene glycol 3,350	10.	7.0	
11.	0.2 M Potassium iodide	11.	20% w/v Polyethylene glycol 3,350	11.	7.0	
12.	0.2 M Ammonium iodide	12.	20% w/v Polyethylene glycol 3,350	12.	6.2	
13.	0.2 M Sodium thiocyanate	13.	20% w/v Polyethylene glycol 3,350	13.	6.9	
14.	0.2 M Potassium thiocyanate	14.	20% w/v Polyethylene glycol 3,350	14.	7.0	
15.	0.2 M Lithium nitrate	15.	20% w/v Polyethylene glycol 3,350	15.	7.1	
16.	0.2 M Magnesium nitrate hexahydrate	16.	20% w/v Polyethylene glycol 3,350	16.	5.9	
17.	0.2 M Sodium nitrate	17.	20% w/v Polyethylene glycol 3,350	17.	6.8	
18.	0.2 M Potassium nitrate	18.	20% w/v Polyethylene glycol 3,350	18.	6.8	
19.	0.2 M Ammonium nitrate	19.	20% w/v Polyethylene glycol 3,350	19.	6.2	
20.	0.2 M Magnesium formate dihydrate	20.	20% w/v Polyethylene glycol 3,350	20.	7.0	
21.	0.2 M Sodium formate	21.	20% w/v Polyethylene glycol 3,350	21.	7.2	
22.	0.2 M Potassium formate	22.	20% w/v Polyethylene glycol 3,350	22.	7.3	
23.	0.2 M Ammonium formate	23.	20% w/v Polyethylene glycol 3,350	23.	6.6	
24.	0.2 M Lithium acetate dihydrate	24.	20% w/v Polyethylene glycol 3,350	24.	7.9	
25.	0.2 M Magnesium acetate tetrahydrate	25.	20% w/v Polyethylene glycol 3,350	25.	7.9	
26.	0.2 M Zinc acetate dihydrate	26.	20% w/v Polyethylene glycol 3,350	26.	6.4	
27.	0.2 M Sodium acetate trihydrate	27.	20% w/v Polyethylene glycol 3,350	27.	8.0	
28.	0.2 M Calcium acetate hydrate	28.	20% w/v Polyethylene glycol 3,350	28.	7.5	
29.	0.2 M Potassium acetate	29.	20% w/v Polyethylene glycol 3,350	29.	8.1	
30.	0.2 M Ammonium acetate	30.	20% w/v Polyethylene glycol 3,350	30.	7.1	
31.	0.2 M Lithium sulfate monohydrate	31.	20% w/v Polyethylene glycol 3,350	31.	6.0	
32.	0.2 M Magnesium sulfate heptahydrate	32.	20% w/v Polyethylene glycol 3,350	32.	6.0	
33.	0.2 M Sodium sulfate decahydrate	33.	20% w/v Polyethylene glycol 3,350	33.	6.7	
34.	0.2 M Potassium sulfate	34.	20% w/v Polyethylene glycol 3,350	34.	6.8	
35.	0.2 M Ammonium sulfate	35.	20% w/v Polyethylene glycol 3,350	35.	6.0	
36.	0.2 M Sodium tartrate dibasic dihydrate	36.	20% w/v Polyethylene glycol 3,350	36.	7.3	
37.	0.2 M Potassium sodium tartrate tetrahydrate	37.	20% w/v Polyethylene glycol 3,350	37.	7.4	
38.	0.2 M Ammonium tartrate dibasic	38.	20% w/v Polyethylene glycol 3,350	38.	6.6	
39.	0.2 M Sodium phosphate monobasic monohydrate	39.	20% w/v Polyethylene glycol 3,350	39.	4.7	
40.	0.2 M Sodium phosphate dibasic dihydrate	40.	20% w/v Polyethylene glycol 3,350	40.	9.1	
41.	0.2 M Potassium phosphate monobasic	41.	20% w/v Polyethylene glycol 3,350	41.	4.8	
42.	0.2 M Potassium phosphate dibasic	42.	20% w/v Polyethylene glycol 3,350	42.	9.2	
43.	0.2 M Ammonium phosphate monobasic	43.	20% w/v Polyethylene glycol 3,350	43.	4.6	
44.	0.2 M Ammonium phosphate dibasic	44.	20% w/v Polyethylene glycol 3,350	44.	8.0	
45.	0.2 M Lithium citrate tribasic tetrahydrate	45.	20% w/v Polyethylene glycol 3,350	45.	8.4	
46.	0.2 M Sodium citrate tribasic dihydrate	46.	20% w/v Polyethylene glycol 3,350	46.	8.3	
47.	0.2 M Potassium citrate tribasic monohydrate	47.	20% w/v Polyethylene glycol 3,350	47.	8.3	
48.	0.2 M Ammonium citrate dibasic	48.	20% w/v Polyethylene glycol 3,350	48.	5.1	



\diamond Measured pH at 25 ° C

PEG/Ion Screen contains forty-eight unique reagents. To determine the formulation of each reagent, simply read across the page.

Tube #	Salt	Tube #	Buffer \diamond	Tube #	Polymer
1.	0.1 M Sodium malonate pH 4.0	1.	None	1.	12% w/v Polyethylene glycol 3,350
2.	0.2 M Sodium malonate pH 4.0	2.	None	2.	20% w/v Polyethylene glycol 3,350
3.	0.1 M Sodium malonate pH 5.0	3.	None	3.	12% w/v Polyethylene glycol 3,350
4.	0.2 M Sodium malonate pH 5.0	4.	None	4.	20% w/v Polyethylene glycol 3,350
5.	0.1 M Sodium malonate pH 6.0	5.	None	5.	12% w/v Polyethylene glycol 3,350
6.	0.2 M Sodium malonate pH 6.0	6.	None	6.	20% w/v Polyethylene glycol 3,350
7.	0.1 M Sodium malonate pH 7.0	7.	None	7.	12% w/v Polyethylene glycol 3,350
8.	0.2 M Sodium malonate pH 7.0	8.	None	8.	20% w/v Polyethylene glycol 3,350
9.	4% v/v Tacsimate pH 4.0	9.	None	9.	12% w/v Polyethylene glycol 3,350
10.	8% v/v Tacsimate pH 4.0	10.	None	10.	20% w/v Polyethylene glycol 3,350
11.	4% v/v Tacsimate pH 5.0	11.	None	11.	12% w/v Polyethylene glycol 3,350
12.	8% v/v Tacsimate pH 5.0	12.	None	12.	20% w/v Polyethylene glycol 3,350
13.	4% v/v Tacsimate pH 6.0	13.	None	13.	12% w/v Polyethylene glycol 3,350
14.	8% v/v Tacsimate pH 6.0	14.	None	14.	20% w/v Polyethylene glycol 3,350
15.	4% v/v Tacsimate pH 7.0	15.	None	15.	12% w/v Polyethylene glycol 3,350
16.	8% v/v Tacsimate pH 7.0	16.	None	16.	20% w/v Polyethylene glycol 3,350
17.	4% v/v Tacsimate pH 8.0	17.	None	17.	12% w/v Polyethylene glycol 3,350
18.	8% v/v Tacsimate pH 8.0	18.	None	18.	20% w/v Polyethylene glycol 3,350
19.	0.1 M Succinic acid pH 7.0	19.	None	19.	12% w/v Polyethylene glycol 3,350
20.	0.2 M Succinic acid pH 7.0	20.	None	20.	20% w/v Polyethylene glycol 3,350
21.	0.1 M Ammonium citrate tribasic pH 7.0	21.	None	21.	12% w/v Polyethylene glycol 3,350
22.	0.2 M Ammonium citrate tribasic pH 7.0	22.	None	22.	20% w/v Polyethylene glycol 3,350
23.	0.1 M DL-Malic acid pH 7.0	23.	None	23.	12% w/v Polyethylene glycol 3,350
24.	0.2 M DL-Malic acid pH 7.0	24.	None	24.	20% w/v Polyethylene glycol 3,350
25.	0.1 M Sodium acetate trihydrate pH 7.0	25.	None	25.	12% w/v Polyethylene glycol 3,350
26.	0.2 M Sodium acetate trihydrate pH 7.0	26.	None	26.	20% w/v Polyethylene glycol 3,350
27.	0.1 M Sodium formate pH 7.0	27.	None	27.	12% w/v Polyethylene glycol 3,350
28.	0.2 M Sodium formate pH 7.0	28.	None	28.	20% w/v Polyethylene glycol 3,350
29.	0.1 M Ammonium tartrate dibasic pH 7.0	29.	None	29.	12% w/v Polyethylene glycol 3,350
30.	0.2 M Ammonium tartrate dibasic pH 7.0	30.	None	30.	20% w/v Polyethylene glycol 3,350
31.	2% v/v Tacsimate pH 4.0	31.	0.1 M Sodium acetate trihydrate pH 4.6	31.	16% w/v Polyethylene glycol 3,350
32.	2% v/v Tacsimate pH 5.0	32.	0.1 M Sodium citrate tribasic dihydrate pH 5.6	32.	16% w/v Polyethylene glycol 3,350
33.	2% v/v Tacsimate pH 6.0	33.	0.1 M BIS-TRIS pH 6.5	33.	20% w/v Polyethylene glycol 3,350
34.	2% v/v Tacsimate pH 7.0	34.	0.1 M HEPES pH 7.5	34.	20% w/v Polyethylene glycol 3,350
35.	2% v/v Tacsimate pH 8.0	35.	0.1 M Tris pH 8.5	35.	16% w/v Polyethylene glycol 3,350
36.	None	36.	0.07 M Citric acid, 0.03 M BIS-TRIS propane / pH 3.4	36.	16% w/v Polyethylene glycol 3,350
37.	None	37.	0.06 M Citric acid, 0.04 M BIS-TRIS propane / pH 4.1	37.	16% w/v Polyethylene glycol 3,350
38.	None	38.	0.05 M Citric acid, 0.05 M BIS-TRIS propane / pH 5.0	38.	16% w/v Polyethylene glycol 3,350
39.	None	39.	0.04 M Citric acid, 0.06 M BIS-TRIS propane / pH 6.4	39.	20% w/v Polyethylene glycol 3,350
40.	None	40.	0.03 M Citric acid, 0.07 M BIS-TRIS propane / pH 7.6	40.	20% w/v Polyethylene glycol 3,350
41.	None	41.	0.02 M Citric acid, 0.08 M BIS-TRIS propane / pH 8.8	41.	16% w/v Polyethylene glycol 3,350
42.	0.02 M Calcium chloride dihydrate, 0.02 M Cadmium chloride hydrate, 0.02 M Cobalt(II) chloride hexahydrate	42.	None	42.	20% w/v Polyethylene glycol 3,350
43.	0.01 M Magnesium chloride hexahydrate 0.005 M Nickel(II) chloride hexahydrate	43.	0.1 M HEPES sodium pH 7.0	43.	15% w/v Polyethylene glycol 3,350
44.	0.02 M Zinc chloride	44.	None	44.	20% w/v Polyethylene glycol 3,350
45.	0.15 M Cesium chloride	45.	None	45.	15% w/v Polyethylene glycol 3,350
46.	0.2 M Sodium bromide	46.	None	46.	20% w/v Polyethylene glycol 3,350
47.	1% w/v Tryptone	47.	0.05 M HEPES sodium pH 7.0	47.	12% w/v Polyethylene glycol 3,350
48.	1% w/v Tryptone	48.	0.05 M HEPES sodium pH 7.0	48.	20% w/v Polyethylene glycol 3,350

\diamond Buffer pH is that of a 1.0 M stock prior to dilution
with other reagent components: pH with HCl or NaOH.

PEG/Ion 2 Screen contains forty-eight unique reagents. To determine the formulation of each reagent, simply read across the page.

34 Journey

Aliso Viejo, CA 92656-3317 U.S.A.

Tel: (949) 425-1321 • Fax: (949) 425-1611

E-mail: tech@hrmail.com

Website: www.hamptonresearch.com

HAMPTON
RESEARCH

Solutions for Crystal Growth

© 1991-2010 Hampton Research Corp. all rights reserved
Printed in the United States of America. This guide or
parts thereof may not be reproduced in any form without
the written permission of the publishers.