

Miltenyi MacsQuant VRB, 3 lasers/8 fluorescences + passeur de plaques, L2, Bâtiment CIRI-SFR, Gerland 4B-2R-2V



Laser	Channel	Filter (nm)	Dye or parameter
Violet 405 nm	V1	450/50	CFP, VioBlue
	V2	525/50	Pacific Orange™, VioGreen
Blue 488 nm	B1	525/50	GFP, FITC
	B2	585/40	PE
	B3	655-730	PI, PerCP, PE-Cy™5.5, PerCP-Vio700, PE-Vio615, ECD, PE-CF594, PE/Dazzle™ 594, PE-eFluor® 610
	B4	750 LP	PE-Cy7, PE-Vio770
Red 635 nm	R1	655-730	APC
	R2	750 LP	APC-Cy7, APC-Vio770
Blue 488 nm	FSC	488/10	Size
	SSC	488/10	Granularity